

# Material Safety Data Sheet

Identity	Scrip Ultrasound Lotion
----------	-------------------------

## HAZARDOUS INGREDIENTS

Ingredients	%	CAS Number	LD/50 (Species & Route)	LC/80 (Species & Route)
<b>THIS PRODUCT IS KNOWN TO BE A NON-HAZARDOUS PRODUCT</b>				

## PHYSICAL DATA

Physical State	Liquid	Vapour Pressure	n/a
Specific Gravity (Density)	0.97	Vapour Density	n/a
Boiling Point (°C)	100	Percent Volatile by Volume	n/a
Freezing Point (°C)	0	Evaporation Rate	n/a
pH	6.50-7.50	Solubility in Water (20°C)	good
Odour & Appearance	Slight fragrance. White to creamy white emulsion.		
Odour Threshold	n/a		

## FIRE AND EXPLOSION HAZARD DATA

Flammability Classification	Non-flammable
Flammable Limits	
Lower Explosion Level (% by Volume)	n/a
Upper Explosion Level (% by Volume)	n/a
Flashpoint	n/a
Extinguish With	n/a

# Material Safety Data Sheet

## REACTIVITY DATA

Chemical Stability	n/a
Incompatible Materials	n/a
Hazardous Polymerization	n/a
Hazardous Decomposition	n/a

## TOXICOLOGICAL PROPERTIES

Route of Entry	Symptoms
Skin Contact	n/a
Skin Absorption	n/a
Eye Contact	May cause redness to eye.
Inhalation Acute	n/a
Inhalation Chronic	n/a
Ingestion	May cause stomach upset if ingested.
Effects of Chronic Exposure to Product	n/a
Irritancy of Product	n/a
Exposure Limits of Product	n/a
Sunsilization to Product	n/a
Synergistic Materials	n/a

## PREVENTATIVE MEASURES

Engineering Controls	n/a
Waste Disposal	Dispose of in accordance with local, state and federal regulations.
Leak and Spill Procedure	Scoop up excess. Clean with appropriate cleaner and rinse with water.
Handling Procedures & Equipment	No special requirements.
Storage Requirements	Store away from heat and extreme cold.
Special Shipping Instructions	Should not freeze.

# Material Safety Data Sheet

---

## FIRST AID MEASURES

Skin Contact	n/a
Eye Contact	Flush with water for 15 minutes while holding eyelids open. Seek medical attention if irritation persists.
Inhalation	n/a
Ingestion	Induce vomiting. Seek medical attention.