

SAFETY DATA SHEET

Issuing Date 12-Dec-2013

Revision Date 10-Jun-2015

Revision Number 1

1. IDENTIFICATION OF THE SUBSTANCE/PREPARATION AND THE COMPANY/UNDERTAKING

GHS product identifier

Product Name Cardinal Health Instant Cold Packs, Blue Color

Other means of identification

UN-Number UN1942

Synonyms Cardinal Health Porta Sample Cold Pack
Kwik-Kold™ Instant Cold Packs

Recommended use of the chemical and restrictions on use

Recommended Use Instant Cold Pack for thermal therapy

Uses advised against No information available

Supplier's details

Supplier Address

Cardinal Health
7000 Cardinal Place
Dublin, OH 43017
TEL: 1-800-964-5227 (Customer Service, 7:30 am (EST) – 8:00 pm (EST))

Emergency telephone number

Emergency Telephone Number For emergency information Contact: CHEMTREC
For the United States and Canada: 1-800-424-9300
For international shipments: +1 703-527-3887

2. HAZARDS IDENTIFICATION

Classification

If container is damaged will release a chemical considered hazardous according to the OSHA Hazard Communication Standard 2012 (29 CFR 1910.1200)

A waterproof package containing small white ammonium nitrate pellets surrounding a small rupturable plastic bag filled with water. Once contents are mixed, the bag will feel cold to the touch. Used as a first aid device. No hazards expected when used per instructions. The hazards described below apply if the packaging is damaged resulting in release of ammonium nitrate solution.

Serious Eye Damage/Eye Irritation	Category 2A
Oxidizing solids	Category 3

GHS Label elements, including precautionary statements

Emergency Overview

Signal Word Hazard Statements <ul style="list-style-type: none">• Causes serious eye irritation• May intensify fire; oxidizer	Warning
---	----------------