

EASY TAP™

Urinary Leg Bags

Your Solution for Easy, Simple & Effective Urinary Containment.

Comfortable Soft Leg Bag with an Easy to Use Drain Valve

The Rüsche Easy Tap™ Urinary Leg Bag provides a simple and comfortable solution. The Easy Tap™ is a high-quality vinyl bag with an anti-reflux valve that reduces the risk of backflow, which reduces the danger of bacteria and infection. It is manufactured with leak proof, electronically welded seams.

The Easy Tap™ Flip Valve is an improved drain valve enabling convenient one-handed operation even with limited dexterity. Provides easy and safe drainage.

Anti-Reflux Valve

Reduces the risk of backflow, which reduces the danger of infection. Has a longer catheter adapter for greater security with all types of catheters. Fits all latex & PVC extension tubing with its universal connector.

High-Quality Vinyl Bag

Comfortable and quiet. Soft vinyl bag with rounded edges and easy to read volume markings. Provides fluid sterile pathway. Can be worn on the inside thigh or calf.

EasyTap™ Flip Valve

Improved drain valve design that enables convenient, one-handed operation. The valve is both easy and smooth to open and close. There is no bottom cap to lose or misplace.




Adjustable fabric straps

Soft, breathable material reducing the risk of irritation to the leg. Provides a comfortable method for securing the urinary leg bag to the user, eliminating the problem of slippage and disconnection. Made with reinforced notches and buttons for easy attachment. Can be easily cut to exact size needed.

EASY TAP™ LEG BAGS 48/Box

Product #	Volume	Size
452919	500 ML	Medium
452932	1000 ML	Large



To place an order contact
Customer Service at **1-800-553-5214**

2917 Weck Drive
Research Triangle Park, NC • 27709
www.rusch.com
P: 800 553 5214 • F: 770 623 1829

RÜSCH
A TELEFLEX COMPANY