

SuperEar®

SE5000

Personal Sound Amplifier

HANDHELD SOUND AMPLIFIER

Automatic Shut Off Microphone	No
Body Frame Material	50dB gain, multielement, swivel type
Max Output	High Impact ABS
Humidity Operating Range	107 dB
Operating Temperature	At 20°C, 95% RH
Frequency Range	0°C - 55°C
Plug-in Connector	100Hz - 10kHz
Power Source	Stereo Jack, 3.5mm diameter
Weight	One AAA Alkaline battery (included)
Belt/Pocket Clip	3 oz.
	Included

HEADPHONES (SE-HP)

Diaphragm	Mylar, diameter 30mm.
Impedance	32 ohms +/- 15% at 1kHz
Nominal Power	4mW
Max Power	10mW
Sensitivity	100 dB S.P.L +/- 2dB 1kHz
Balance per Channel	+/- 3 dB at 1kHz
Frequency Response	20Hz-20kHz
Plug-in Connector	Stereo Plug, 3.5mm diameter
Cord Length	1200mm - +/-50mm
Weight	50g +/-5g
Earpad	Black, cushioned sponge

EARBUDS (LC-EB)

Diaphragm	Mylar, diameter 10mm.
Impedance	32 ohms +/- 15% at 1kHz
Nominal Power	5mW
Max Power	10mW
Sensitivity	105 dB S.P.L +/- 3dB 1kHz
Balance per Channel	+/- 3 dB at 1kHz
Frequency Response	20Hz-20kHz 1kHz
Plug-in Connector	Stereo Plug, 3.5mm diameter
Cord Length	1150mm - +/-50mm
THD	Less than 5% at 1kHz 0.126V
Weight	15g +/-5g



HS TARIFF CODE: 8515.50

SE5000 DETAILS

WARRANTY	3 year
MEASUREMENTS	3.05"x 2"x 0.75"
CERTIFIED	CE & ROHS
COMPLIANT	CMS, MDS, ADA, ACA 1557, Prop 65

Cleaning instructions available



SonicTechnology®

Quality and Innovation since 1981

Use with Disposable Sanitary Earpad Covers
- Model EPC100

www.sonictechnology.com
superear@sonictechnology.com

FAX 530-272-4257 | 1-800-247-5548
P.O. Box 539, Grass Valley, CA 95945

V1.11/11/21/19