

2320012 / 2320185

Material Safety Data Sheet
For
CryoStim

I. PRODUCT IDENTIFICATION

TRADE NAME (As Labeled): CryoStim Conductive Gel

PRODUCT CLASS AND/OR

PRODUCT DESIGNATION: A LUBRICANT GEL TO BE USED WITH THE CRYOSTIM PROBE.

DATE PREPARED OR REVISED: MAY 10, 2002

II. HAZARDOUS INGREDIENTS

NO INGREDIENTS KNOWN TO BE HAZARDOUS ARE PRESENT. THE IDENTITY OF THE INDIVIDUAL COMPONENTS IS PROPRIETARY INFORMATION AND IS REGARDED TO BE A TRADE SECRET. IN ACCORDANCE WITH 29 CFR 1910.1200 WE HAVE TO ASSUME THAT THE MIXTURE PRESENTS THE SAME HEALTH HAZARDS AS THE INDIVIDUAL COMPONENTS WHEN THEY ARE PRESENT AT GREATER THAN 10% CONCENTRATION.

III. PHYSICAL PROPERTIES

FLASH POINT: NOT DETERMINED
BOILING POINT: NOT DETERMINED
EVAPORATION: NOT DETERMINED
ETHER
VAPOR DENSITY: NOT DETERMINED

VAPOR PRESSURE: NOT DETERMINED
pH: (mmHg) 5.8
DENSITY: 8.0 (lbs/gal)
SOLUBILITY IN WATER: SOLUBLE
APPEARANCE: A LIGHT GREEN GEL WITH AN HERBAL FRAGRANCE.

IV. FIRE AND EXPLOSION

EXTINGUISHING MEDIA: A NON-FLAMMABLE LIQUID GEL.
UNUSUAL FIRE AND EXPLOSION HAZARDS: NONE
SPECIAL FIRE FIGHTING PROCEDURES: NONE REQUIRED
STABILITY: STABLE
INCOMPATIBILITY: NOT DETERMINED
HAZARDOUS DECOMPOSITION PRODUCTS: SODIUM CHLORIDE