

232 0214

MATERIAL SAFETY DATA SHEET

I. IDENTITY:

Label Name: **Cramer Flex-I-Cold Packs**
Item Numbers: 032746, 032846

Date Prepared: 11/1/06

Chemical Name and Synonyms: N A
Chemical Family: Mixture

Package types: Flexible pouch
Product Class: Reuseable Cold Pack

II. HAZARDOUS INGREDIENTS:

Hazardous Components	CAS #	OSHA	PEL	ACGIH	TLV	%
None						

III. PHYSICAL CHARACTERISTICS:

Boiling Point: Unknown	Specific Gravity (Water=1): ~1
Vapor Pressure: Unknown	Melting Point: ~ -4° F
Vapor Density: (Air=1): Unknown	Evaporation Rate: Unknown
Solubility in Water: Somewhat soluble	Percent Volatile by Volume: ~ 75%
Appearance and Odor: Clear blue, viscous gel	

IV. FIRE AND EXPLOSION HAZARD:

Flash Point: Unknown	Flammable Limits: Unknown
Extinguishing Media: Carbon dioxide, dry chemical or foam	
Special Fire Fighting Procedures: None Known.	
Unusual Fire and Explosion Hazards: None Known.	

V. REACTIVITY DATA:

Stability: Stable	Conditions to Avoid: None known
Incompatibility to Avoid: None known	
Hazardous Decomposition Byproducts: None known	
Hazardous Polymerization: Will not occur at room temperatures	

VI. HEALTH HAZARD DATA:

Eyes: Gel may cause minor temporary eye irritation if it comes into contact with eyes.

Skin: Frostbite or cold injury may occur if chilled pack is applied directly to skin. Follow directions on pack.

Ingestion: Do not take internally. May cause gastrointestinal discomfort and diarrhea.

Inhalation: N/A.

Emergency First Aid Procedures:

Skin: Treat for frostbite or cold injury as appropriate.

Eyes: Flush eyes with water for at least 15 minutes. If irritation persists, or if material was hot and eyes are burned, contact physician.

Ingestion: Induce vomiting. Watch for symptoms of salt (sodium chloride) poisoning.

VII. PRECAUTIONS FOR SAFE HANDLING AND USE:

Spills: Spills should be contained, collected into absorbent material and disposed of in accordance with all applicable federal, state and local regulations.

Handling and Storage: Store at moderate room temperature (50-110 ° F). Read label carefully and follow directions.

VIII. CONTROL MEASURES:

Respiratory Protection: None required.

Ventilation: No special requirements.

Eye Protection: Not required with normal use.

Protective Gloves: Not required.

Protective Clothing or Equipment: Not required.

Abbreviation Key: N/A = not applicable.

The information contained herein is believed to be accurate. It is the user's obligation to determine the safe use of the product.