

## Duet Rollator/Transport Chair, 8" Casters

HCPCS: E0143+E0156




FIGURE A



FIGURE B




FIGURE C




FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Combines the features of a rollator and transport chair in one unit An individual can ambulate independently or be pushed safely by a caregiver
- To use as a transport chair the Duet has locking flip-up armrests (Figure A - Figure C) and a flip-down footrest (Figure D) to aid the user A safety lock release is used to store the arms and a hook and loop fastener strap is used to store the footrest while using the Duet for ambulation
- The tool-free two-position contoured padded backrest is reversible depending on whether you are using the Duet as a rollator or in transport chair mode
- Comes with padded seat with carry pouch for added privacy and security of personal items Inside the pouch is a handle to aid with the folding
- The folding frame has a lock-out feature to ensure the unit stays open while in use and ensures the unit will stay folded during transport
- 8" casters are ideal for indoor and outdoor use
- Handles are adjustable in height, with self-threading knob (Figure E)

### PRODUCTS

#### Description

|       |  |
|-------|--|
| 795B  | Duet Rollator/TransportChair,Blue, 1/ea  |
| 795BK | Duet Rollator/TransportChair,B           |
| 795BU | Duet Rollator/TransportChair,Burgdy,1/cs |

### SPECIFICATIONS

#### Description

|                        |                               |
|------------------------|-------------------------------|
| Back Height            | 35"                           |
| Carton Shipping Weight | 25 lbs.                       |
| Handle Height          | 31.5" - 37"                   |
| Overall Dimensions     | 25.5" (W) x 37" (H) x 26" (L) |
| Seat Dimensions        | 13.5" (W) x 13" (D)           |
| Seat To Floor Height   | 21"                           |
| Seat Width             | 13.5"                         |
| Warranty               | Limited Lifetime              |
| Weight                 | 20 lbs.                       |
| Weight Capacity        | 300 lbs                       |

### PARTS

#### Description

|            |  |
|------------|--|
| 9501A79512 | Handlebar, Right, 795, 1/ea, EA            |
| 9502M7951L | Elbow Rest, Left, 795, 2M Ser.#, 1/ea, EA  |
| 9502M7951R | Elbow Rest, Right, 795, 2M Ser.#, 1/ea, EA |
| 750FORK-2M | Front Fork, 795, 2M Ser.#, 1/ea, EA        |

| PARTS         | Description                                  |
|---------------|--|
| 9501A7953BK   | FootrestAssm.Black,795 1A/2Y Ser.#,1ea, EA   |
| 9501A7953R-2M | FootrestAssembly,Red,795 2M Ser.#, 1/ea, EA  |
| 9501A79511    | Handlebar, Left, 795, 1/ea, EA               |
| 9501A7953B-2M | FootrestAssembly,Blue,795 2M Ser.#, 1/ea, EA |
| 9501A79510-RP | Roll Pin Fork Fastener,795,1ea, EA           |
| 9501A79510-2M | Fork Stem Bolt, 795, 2M Ser.#, 1/ea, EA      |
| 9501A79506    | RearAxleBolt795,1ea, EA                      |
| 9501A79506-2M | Rear Axle Bolt, 795, 2M Ser.#, 1/ea, EA      |
| 9501A79504    | Tote for 795 Rollator, 1/ea, EA              |
| 750BCL-2M     | HandBrakeAssmby,Left,795,2M Ser.#, 1/ea, EA  |
| 750BRAKE-2M   | Brake, Right or Left, 795, 2M S#, 1ea, EA    |
| 9501A75003    | Spring Boot for Cable,1/ea, EA               |
| 9500075005    | Folding Hinge Bolt,750,1/ea, EA              |
| 9501A75013    | Seat Bracket, 795 Rollator, 2/pr, EA         |
| 9502M75013    | SeatBrackets&Screws,795 2M Ser.#, 1/pr, PR   |
| 9501A75008    | 3 Hole Wheel,8"Rear,750N,795,1/ea, EA        |