



wenzelite



HCPCS  
E1236  
E1226

Accommodates users with limited hip flexion, allows for gravity-assisted positioning for improved visual orientation, speech and swallowing and aids with bowel and bladder management.

## Pediatric Viper Plus Reclining Wheelchair

Deluxe High Strength, Lightweight

Item #	Size	Color	Overall Width (open)	Overall Width (folded)	Seat Dimensions	Product Weight	Carton Shipping Weight	UOM
PL412RBDDA	12"	Black	20"	13"	12" W x 12" D	42 lbs.	48 lbs.	1/bx
PL414RBDDA	14"	Black	22.5"	13.5"	14" W x 14" D	43 lbs.	49 lbs.	1/bx

- State-of-the art hydraulic reclining mechanism.
- Backrest reclines from 90° to 180°.
- Double-coated chrome carbon steel frame for maximum durability.
- Durable flame retardant nylon upholstery.
- Padded flip-back detachable desk armrests.
- Composite lightweight 24" Mag-style rear wheels with urethane tires and 6" front rubber tires.
- Offset rear wheel position provides additional stability and evenly distributes user's weight in all reclining positions.
- Precision sealed front and rear wheel bearings.
- Front caster forks are adjustable in three positions.
- Standard dual axle provides easy transition of seat to hemi-level.
- Push-to-lock wheel locks.
- Spreader bar functions as a back stabilizer and a push handle for caregiver.
- Standard with height adjustable swing-away elevating leg rests.
- Durable and lightweight composite footplates.
- Flip-up anti-tippers with wheels.

### specifications

Seat-To-Floor Height	15", 16"
Back Height	26"
Overall Length w/Riggings	38"
Overall Height	40.5"
Weight Capacity	250 lbs.
Warranty	Limited Lifetime

### accessories

STDS850	Seat belt, auto clasp type	1/bx
STDS851	Seat belt, Velcro type closure	1/bx
STDS810	2" seat rail extension	1pr/bx